


THIS SET OF FINAL DRAWINGS IS REFLECTIVE OF THE LATEST SALES CONFIGURATION. ANY CHANGES TO THIS SALES CONFIGURATION MAY REQUIRE A REVISION TO THIS PROJECT PLAN. IF REQUESTED, SIEMENS WILL PRODUCE A REVISED SET OF FINAL DRAWINGS TO REFLECT THE CHANGES, HOWEVER SIEMENS IS NOT RESPONSIBLE FOR ANY CONSTRUCTION COSTS ASSOCIATED WITH THE CHANGES THAT OCCUR FROM THIS PLAN MODIFICATION.





STRETCHER SIZE (6'-5" X 2'-7") SHOWN IS FOR REFERENCE ONLY. VERIFICATION AND COORDINATION BY CUSTOMER IS REQUIRED TO ENSURE PROPER TRANSPORT AND WORKFLOW ACCESS.

HEIGHT OF WINDOW TO BE COORDINATED WITH COUNTERTOP/DESK HEIGHT.

### ARCHITECTURAL EQUIPMENT PLAN

SCALE: 1/4" = 1'-0"



**STATE AGENCY REVIEW**

PRIOR TO SIEMENS EQUIPMENT INSTALLATION, APPROVAL OF CONSTRUCTION OR STRUCTURAL MODIFICATIONS UTILIZING X-RAY FOR DIAGNOSTIC OR THERAPEUTIC PURPOSES, MUST BE OBTAINED BY THE CUSTOMER FROM THE APPROPRIATE STATE AGENCY, IF APPLICABLE.

FINISHED ROOM HEIGHT	
FOR CT GANTRY ONLY	MINIMUM 6'-10 11/16"
FOR CT GANTRY WITH INJECTOR ARM	MINIMUM 7'-6 9/16"
CAREVISION MONITOR/CEILING MOUNT	SEE DETAIL ON S-102 SHEET

SYM	DATE	DESCRIPTION
02/04/26	R-101RA VERSION DATED 02/04/26 APPROVED BY CUSTOMER FOR FINALS	

PROJECT MANAGER: MIKE ROWAN  
 TEL: (813) 431-2205 EXT:  
 FAX: (813) 315-6350  
 EMAIL: MIKE.ROWAN@SIEMENS-HEALTHINEERS.COM

**SIEMENS**

**ORLANDO HEALTH / MELBOURNE FSED**  
 LAKE ANDREW DRIVE, MELBOURNE, FL 32940  
 CT SCAN ROOM - SOMATOM GO-TOP

PROJECT #: **2610370** SHEET: **A-101**

DATE: T. KELLEY  
 DRAWN BY: T. KELLEY

SCALE: AS NOTED REF. # CPQ-1648726/P

Facility Name: Orlando Health Melbourne Hospital  
 Project Name: Freestanding ER - Viera  
 OPC Project #: 2323980034-16987

### ATTENTION:

- THIS DRAWING IS DESIGNED TO CONFORM TO FEATURES AND EQUIPMENT REQUIREMENTS PRESENTED AT THE TIME OF THEIR PREPARATION. SINCE BOTH THESE FACTORS ARE SUBJECT TO DESIGN MODIFICATION, THEY ARE NOT TO BE USED FOR CONSTRUCTION PURPOSES.

- IT IS RECOMMENDED THAT THE SIEMENS DRAWINGS BE INCORPORATED WITH THE CONSTRUCTION DOCUMENTS FOR REFERENCE.

- ALL DIMENSIONS SHOWN ON THIS DRAWING ARE FROM FINISHED SURFACES.  
 - THIS DRAWING DOES NOT PROVIDE RADIATION SHIELDING REQUIREMENTS FOR X-RAY AND ASSOCIATED EQUIPMENT. THE CUSTOMER IS RESPONSIBLE FOR CONSULTING WITH A REGISTERED RADIATION PHYSICIST TO SPECIFY RADIATION PROTECTION.

### EQUIPMENT LEGEND

NO	DESCRIPTION	SMS SYM	WEIGHT (LBS)	BTU/HR TO AIR	DIMENSIONS (INCHES)			REMARKS
					W	D	H	
1	23" FLAT SCREEN WITH DUAL MONITORS, KEYBOARD AND CONTROL DEVICE	Ⓢ	20	---	22 1/2	9 1/4	19 3/8	
2	SOMATOM GO.TOP GANTRY	Ⓢ	2,855	24,226	87 1/8	32 3/4	73 1/8	
3	PATIENT TABLE - VARIO 2	Ⓢ	780	1,024	99 13/16	27 3/8	21 13/16	307kg, MAX. HEIGHT 38"
4	WIRELESS ACCESS POINT	Ⓢ	2.5	---	---	---	---	HEIGHT OFF FF: > 6'-6 3/4"
5	POWER OVER ETHERNET SWITCH	Ⓢ	---	---	6 7/16	2	1 3/16	
6	EATON SURGE PROTECTIVE DEVICE PANEL	Ⓢ	13.5	---	7 1/2	6 11/16	12	WALL MOUNTED

### PROJECT MILESTONES TO BE COMPLETED BEFORE EQUIPMENT DELIVERY

CHECK STATUS	COMPLETION DATE	MILESTONES	REFERENCE SHEET
		SYSTEM STANDARDS	
		ARCHITECTURAL	
		STORAGE AREA AVAILABLE FOR STORING ITEMS DURING INSTALLATION.	A-101/A-102/A-501
		DELIVERY PATH VERIFIED. DELIVERY PATH MUST BE APPROVED BY CUSTOMER / SeOR AND DOCUMENTED.	A-101/A-102/A-501
		MINIMUM DOOR OPENING.	A-101/A-102/A-501
		CLIMATE CONTROL FUNCTIONING 24 HOURS A DAY, 7 DAYS A WEEK.	A-101
		CASEWORK COMPLETE IN EXAM AND CONTROL ROOMS.	A-101
		ROOM LIGHTING COMPLETE AND FUNCTIONING.	A-101
		ALL ROOMS CONTAINING SIEMENS EQUIPMENT ARE CLEAN AND DUST FREE.	A-101
		LEAD SHIELDING (WALLS, DOORS, WINDOWS) COMPLETE.	A-101/A-102
		NETWORK ADDRESSES OBTAINED FOR SIEMENS REMOTE SERVICES (SRS).	A-102/E-102/E-501
		NETWORK DATA DROPS INSTALLED (LANDED AND ACTIVE).	A-102/E-102/E-501
		I.P. AND PACS INFORMATION AVAILABLE AND PROVIDED TO SIEMENS REMOTE SERVICE TEAM (SRS).	A-102/E-102/E-501
		B2B VPN TUNNEL ESTABLISHED BY SRS TEAM.	A-102/E-102/E-501
		NETWORK COORDINATION BETWEEN SITE AND IT SERVICE CONSULTANT COMPLETE.	A-102/E-102/E-501
		STRUCTURAL	
		FLOOR THICKNESS VERIFIED AND WITHIN SPECIFICATIONS.	S-101/S-501
		FOR ANY OEMS, VERIFY IF BUNDLING WITH SYSTEM DELIVERY OR PRE-DELIVERED/DROP-SHIPPED.	S-101/S-102/S-501
		CEILING HEIGHT VERIFIED (CHECK MIN. HEIGHT WITH OPTIONS).	S-102/S-501
		FLOOR LEVELNESS VERIFIED AND WITHIN SPECIFICATIONS.	S-101
		ELECTRICAL	
		CABLES RUNS CHECKED TO ENSURE MAXIMUM LENGTH IS NOT EXCEEDED.	E-101
		ALL CONDUITS, TROUGHS, AND CORE DRILLS ARE OUTSIDE OF THE NO CORE DRILL AREAS.	E-102/E-501
		CABLES INLETS INSTALLED AT LOCATIONS PER PLANS.	E-101/E-102
		MAIN PANEL AND BREAKERS INSTALLED.	E-101/E-102
		CONTRACTOR SUPPLIED ELECTRICAL CABLING AND PIGTAILS INSTALLED.	E-101/E-102
		CONTRACTOR SUPPLIED EPOS INSTALLED, TESTED, AND FUNCTIONING.	E-101/E-102/E-501
		CONTRACTOR SUPPLIED X-RAY WARNING LIGHT AND WIRING INSTALLED.	E-101/E-102/E-501
		CONTRACTOR SUPPLIED DOOR SWITCH AND WIRING INSTALLED (IF APPLICABLE).	E-101/E-102/E-501
		MECHANICAL	
		SYSTEM OPTIONS	
		ACCESSORIES OPTIONS	

### ARCHITECTURAL NOTES

- ALL PRELIMINARY EQUIPMENT LAYOUTS SUBMITTED BY SIEMENS HEALTHCARE ARE BASED ON THE RECOMMENDED SPACE NECESSARY FOR THE OPERATION AND SERVICEABILITY OF THE EQUIPMENT BEING PROPOSED. SIEMENS WILL NOT SUBMIT AN EQUIPMENT LAYOUT THAT IS NOT IN THE BEST INTEREST OF BOTH THE CUSTOMER AND SIEMENS. ALL EQUIPMENT LAYOUTS ARE BASED EITHER ON AN ACTUAL SITE SURVEY OR ARCHITECTURAL DRAWINGS SUPPLIED TO SIEMENS. SIEMENS WILL NOT BE RESPONSIBLE FOR ANY ALTERATIONS THAT ENDOURCH WITHIN DESIGNATED SAFETY AND SERVICE CLEARANCE ZONES AS INDICATED ON DRAWINGS (I.E. PIPE CHASES, VENTILATION DUCTS, CASEWORK, AND SOFFITS, ETC.) MADE BY THE CUSTOMER OR REQUIRED BY A CUSTOMER'S ARCHITECTURAL FIRM ONCE PRELIMINARY DRAWINGS HAVE BEEN SUBMITTED AND APPROVED. DO NOT ALTER ANY SPECIFICATIONS AND/OR DIMENSIONS WITHOUT CONTACTING AND RECEIVING WRITTEN CONFIRMATION FROM SIEMENS PROJECT MANAGER.
- SIEMENS HEALTHCARE IS NOT AN ARCHITECTURAL OR ENGINEERING FIRM. DRAWINGS SUPPLIED BY SIEMENS ARE NOT CONSTRUCTION DRAWINGS. THEREFORE, THESE DRAWINGS ARE TO BE USED ONLY FOR INFORMATION TO COMPLEMENT ACTUAL CONSTRUCTION DRAWINGS AVAILABLE FROM A CUSTOMER APPOINTED ARCHITECTURAL REPRESENTATIVE OR A CUSTOMER'S ENGINEERING DESIGN GROUP. THE CUSTOMER'S ARCHITECT AND GENERAL CONTRACTOR SHALL BE ULTIMATELY RESPONSIBLE FOR COMPLIANCE WITH ALL APPLICABLE CODES AND PROFESSIONAL DESIGN REQUIREMENTS INCLUDING OSHA/NEC SAFETY CLEARANCE REQUIREMENTS IN ADDITION TO SIEMENS-REQUIRED SAFETY/SERVICE CLEARANCES SHOWN.
- THE CUSTOMER IS RESPONSIBLE FOR ALL ROOM AND AREA PREPARATION COSTS, PROFESSIONAL FEES, PERMITS, REPORTS, AND INSPECTION FEES.
- EQUIPMENT WARRANTIES, EXPRESSED OR IMPLIED ON THE PART OF SIEMENS SHALL BE CONTINGENT UPON STRICT COMPLIANCE WITH THE ARCHITECTURAL, STRUCTURAL, ELECTRICAL, MECHANICAL AND RECOMMENDATIONS AND REQUIREMENTS CONTAINED IN THESE DRAWINGS, UNLESS SPECIFIED OTHERWISE.
- ALL DIMENSIONS SHOWN ARE FROM FINISHED SURFACES UNLESS SPECIFIED OTHERWISE.
- THIS DRAWING DOES NOT PROVIDE RADIATION SHIELDING REQUIREMENTS FOR X-RAY AND ASSOCIATED EQUIPMENT. THE CUSTOMER IS RESPONSIBLE FOR CONSULTING WITH A REGISTERED RADIATION PHYSICIST. ACTUAL PROTECTION REQUIREMENTS SHALL BE SPECIFIED BY A REGISTERED RADIATION PHYSICIST AT CUSTOMER'S ENGAGEMENT AND EXPENSE. RESPONSIBILITY FOR ALL INFORMATION AS TO THE ROOM LOCATION, USE, AND NUMBER OF ANTICIPATED EXAMINATIONS TO BE PERFORMED PER TIME PERIOD SHALL BE PROVIDED TO THE PHYSICIST BY THE CUSTOMER. THE CUSTOMER SHALL FURTHER TAKE ALL RESPONSIBILITY IN THE COMMUNICATION AND COORDINATION OF ACTIVITIES OF THE RADIATION PHYSICIST AND THE ARCHITECTURAL REPRESENTATIVE.
- SIEMENS HEALTHCARE SHALL BE RESPONSIBLE FOR SIEMENS EQUIPMENT INSTALLATION, CALIBRATION, CONNECTION AND INSTALLATION OF SIEMENS PROVIDED CABLES. THE CUSTOMER/ELECTRICAL CONTRACTOR IS RESPONSIBLE FOR TERMINATIONS OF CUSTOMER/ELECTRICAL CONTRACTOR-SUPPLIED CABLES TO SIEMENS EQUIPMENT. IN THE EVENT THAT SPECIFIC TRADE RULES OR LICENSE REQUIREMENTS PROHIBIT THIS, THE CUSTOMER SHALL INITIATE THE SERVICES OF APPROVED OTHER CONTRACTORS AND PAY FOR SELECTED, APPROVED PARTIES TO PERFORM THIS WORK WITH SUPERVISION PROVIDED BY SIEMENS. CALIBRATION WHEN ACCOMPLISHED OUTSIDE OF NORMAL INSTALLATION SEQUENCES DUE TO CONTRACTOR OR TRADE RULE ACTIONS OR REQUIREMENTS SHALL BE SUPPORTED BY, CHARGED TO, AND ACCEPTED BY THE CUSTOMER AS AN ADDITIONAL INSTALLATION EXPENSE.
- THE CUSTOMER SHALL COORDINATE WITH SIEMENS PROJECT MANAGER THE LOCATIONS AND TRAVEL OF ALL ANCILLARY EQUIPMENT TO BE CEILING OR WALL MOUNTED (I.E. O.R. LIGHTS, MEDICAL GAS COLUMNS, PHYSIOLOGICAL MONITORING INJECTORS, CRT PLATFORMS, SPRINKLER HEADS, SMOKE DETECTORS, ELECTRICAL OUTLETS, HVAC GRILLES, SPEAKERS, AND GENERAL ROOM LIGHTING, ETC.).
- THE GENERAL CONTRACTOR/CUSTOMER SHALL BE RESPONSIBLE FOR ALL FINAL PAINT, TOUCH-UP AND ANY COSMETIC OR TRIM WORK WHICH NEEDS TO BE OR IS REQUIRED TO BE COMPLETED AFTER THE INSTALLATION OF THE SIEMENS EQUIPMENT AND ANY ASSOCIATED SUPPORT APPARATUS.
- CUSTOMER/CONTRACTOR MUST ASSIST SIEMENS INSTALLERS WITH INSTALLATION OF EQUIPMENT ABOVE 14'-0". REFER TO THE ELECTRICAL NOTES ON SIEMENS SHEET E-101 FOR MORE DETAILS.

### CASEWORK & ACCESSORY NOTES

- ALL CASEWORK IS EITHER EXISTING OR IS TO BE DESIGNED, DETAILED, FURNISHED AND INSTALLED BY THE CUSTOMER AND/OR CONTRACTOR. FOLLOW DESIGN RECOMMENDATIONS INCLUDED HEREWITH, AS THEY ARE ESSENTIAL FOR THE SUCCESSFUL INSTALLATION & OPERATION OF THE SIEMENS EQUIPMENT.
- THE SOUND SYSTEM AND INTERCOM BETWEEN THE EXAMINATION AND CONTROL ROOMS ARE TO BE LOCATED, FURNISHED AND INSTALLED BY THE CUSTOMER/CONTRACTOR.
- ALL FURNITURE (CHAIRS, ETC.) FOR THE CONTROL ROOM ARE TO BE PROVIDED BY THE CUSTOMER.

### PLANNING REQUIREMENTS

EMERGENCY POWER OFF (EPO) BUTTONS ARE REQUIRED IN CONTROL AREA AND AT LEAST ONE LOCATION IN EXAMINATION OR SCAN ROOM.

DOOR (SAFETY) SWITCHES ARE REQUIRED ON ALL DOORS ACCESSING THE EXAMINATION ROOM IN ACCORDANCE WITH LOCAL CODES.

### RESOURCE LIST (SMS USE ONLY)

DESIGNATION	PG NUMBER	DATE
SOMATOM GO	C2-081.891.01.21.02	10.25
COMMON CT	CT00-000.891.04.25.02	07.25
COMMON CT OPTIONS	CT00-000.891.03.58.02	07.25

GO-TOP REV 30

### TRANSPORT AND DELIVERY

#### PACKAGING INFORMATION

COMPONENT	WIDTH	DEPTH	HEIGHT	WEIGHT (WITH TRANSPORT ADAPTERS)	COMPONENT	WIDTH	DEPTH	HEIGHT	WEIGHT
CT GANTRY CRATE	93"	42 1/8"	85 7/16"	3,545 LBS.	CT UPS ONLY	24 5/8"	20 11/16"	9 3/16"	60 LBS.
CT TABLE CRATE	101 3/8"	32 7/16"	59 3/4"	1,305 LBS.					

#### SPECIAL NOTES

- FLOOR LOADING DURING GANTRY TRANSPORT: THE MAXIMUM POSSIBLE FLOOR LOAD (TWO POINT LOAD) PER ROLLER MAY BE REACHED DURING TRANSPORT. IF NECESSARY, COVER THE TRANSPORT ROUTE WITH METAL SHEETS TO DISTRIBUTE THE LOAD.
- TRANSPORT THROUGH A NARROW ENTRANCE: THE GANTRY CAN BE MOVED THROUGH ENTRANCES AS NARROW AS 3'-3 3/8" WHEN THE TRANSPORT ADAPTER IS PARTIALLY REMOVED BUT WITH SUCH A NARROW PASSAGE, SEVERAL COVERS NEED TO BE REMOVED AS WELL. IN SUCH CASES, IT IS ALSO RECOMMENDED TO ARRANGE FOR MORE HELPERS AND OBSERVERS TO MONITOR THE SPACE AROUND THE GANTRY.
- REDUCING THE GANTRY WEIGHT: IN SPECIAL CASES THE GANTRY WEIGHT MUST BE REDUCED BY A FEW POUNDS, E.G. TO USE A FREIGHT ELEVATOR WITH A LIMITED FREIGHT LOAD. THIS CAN BE DONE BY REMOVING CERTAIN PARTS OF THE GANTRY COVERS AS SHOWN IN THE TABLE. NOTE: A SERVICE ENGINEER IS REQUIRED FOR DISASSEMBLY AND REASSEMBLY. STORE THE REMOVED COVERS CAREFULLY AT A SAFE LOCATION TO AVOID DAMAGE DURING TRANSPORT OF THE GANTRY. THE FOLLOWING LIST SHOWS THE WEIGHT OF THE REMOVABLE COVERS: FRONT COVERS: 77 LBS. REAR COVERS: 73 LBS. SIDE COVERS: 31 LBS.

#### TRANSPORT INFORMATION


TOTAL GANTRY TRANSPORT WEIGHT: 3,272 LBS.  
 GANTRY WITHOUT TRANSPORT DEVICE: 2,855 LBS.  
 TRANSPORT DEVICE: 417 LBS.

#### NORMAL TRANSPORT REQUIREMENTS:


DURING THE MOVEMENT OF THE GANTRY THROUGH CORRIDORS THE TRANSPORT CASTERS ARE SWEVELED OUT FOR STABILITY AS SHOWN BELOW. THE MAXIMUM WIDTH IS 5'-4 1/4" AND THE MAX. LENGTH IS 9'-0 3/16" WHEN CASTERS ARE SWEVELED OUT.

#### NARROW SPACE TRANSPORT REQUIREMENTS:

WHEN TRANSPORTING THE GANTRY THROUGH A NARROW SPACE OR DOORWAY, THE TRANSPORT CASTERS ARE SWEVELED IN AS SHOWN BELOW. THE MAXIMUM WIDTH IS 3'-5" AND MAXIMUM LENGTH IS 11'-3 13/16".



#### TRANSPORT DEVICE



### ENVIRONMENTAL CONDITIONS

AIR TEMPERATURE	MINIMUM 64.4°F TO 86°F MAXIMUM
RELATIVE HUMIDITY	20% TO 75%
ABSOLUTE HUMIDITY	MAXIMUM 30 G/M <sup>3</sup> (NO CONDENSATION AT ANY TIME)
TEMPERATURE GRADIENT	MAXIMUM 6 KELVIN PER HOUR
BAROMETRIC PRESSURE	11.6 TO 15.4 PSI
INSTALLATION ALTITUDE	MAXIMUM 6562 FT. A.S.L.

#### AIR QUALITY

EXTERIOR AIR VENTS SHOULD BE EQUIPPED WITH A FILTRATION SYSTEM OF THE FILTER CLASS MERV 8 TO FILTER DUST PARTICLES >10 μm.

THE ROOM AIR SHOULD BE PROTECTED AGAINST CONTAMINATION BY HYDROGEN SULPHIDE, EVEN IN SMALL AMOUNTS. IF A DANGER OF SUCH CONTAMINATION EXISTS, CORRECTIVE ACTIONS HAVE TO BE TAKEN, E.G., EXTRACTOR FANS, SIPHON, MODIFICATION OF VENTILATION INTAKE, ETC.

NO SUNSHINE DIRECTLY ON GANTRY, INSULATION HAS TO BE APPLIED TO WINDOW (EX. CURTAIN)

#### OZONE STERILIZERS

OZONE STERILIZERS SHOULD NOT BE PLANNED OR OPERATED WITH A CT SCANNER IN THE SAME ROOM. IF THE CUSTOMER USES OZONE STERILIZERS IN THE SCANNER ROOM, IT MAY REDUCE THE LIFETIME OF VARIOUS CT COMPONENTS.

Facility Name: Orlando Health Melbourne Hospital  
 Project Name: Freestanding ER - Viera  
 OPC Project #: 2323980034-1987

## ATTENTION:

- THIS DRAWING IS DESIGNED TO CONFORM TO FEATURES AND EQUIPMENT REQUIREMENTS PRESENTED AT THE TIME OF THEIR PREPARATION. SINCE BOTH THESE FACTORS ARE SUBJECT TO DESIGN MODIFICATION, THEY ARE NOT TO BE USED FOR CONSTRUCTION PURPOSES.  
 - THIS SET OF PLANS REPRESENTS A COMPLETE SET OF DETAILS AND SHOULD NOT BE SEPARATED.

- IT IS RECOMMENDED THAT THE SIEMENS DRAWINGS BE INCORPORATED WITH THE CONSTRUCTION DOCUMENTS FOR REFERENCE.


- ALL DIMENSIONS SHOWN ON THIS DRAWING ARE FROM FINISHED SURFACES.  
 - THIS DRAWING DOES NOT PROVIDE RADIATION SHIELDING REQUIREMENTS FOR X-RAY AND ASSOCIATED EQUIPMENT. THE CUSTOMER IS RESPONSIBLE FOR CONSULTING WITH A REGISTERED RADIATION PHYSICIST TO SPECIFY RADIATION PROTECTION.

GO-TOP  
REV 30

		PROJECT MANAGER: MIKE ROWAN TEL: (813) 431-2205 EXT: FAX: (813) 315-6350 EMAIL: MIKE.ROWAN@SIEMENS-HEALTHINEERS.COM		<b>SIEMENS</b>	
		<b>ORLANDO HEALTH / MELBOURNE FSED</b> LAKE ANDREW DRIVE, MELBOURNE, FL 32940 CT SCAN ROOM - SOMATOM GO-TOP			
		PROJECT #: <b>2610370</b>		SHEET: <b>A-501</b>	
		SHEET 2 OF 6		DRAWN BY: T. KELLEY	
		ALL RIGHTS ARE RESERVED.			
		SCALE: AS NOTED	REF. #: CPQ-1648726/P	DATE: T. KELLEY	
		-ISSUE BLOCK-			




SYM	DATE	DESCRIPTION
△	02/04/26	R-1019A VERSION DATED 02/04/26 APPROVED BY CUSTOMER FOR FINALS

NOTE: FOR THE WEIGHTS OF ALL SIEMENS EQUIPMENT SHOWN ON THIS PLAN, SEE THE "EQUIPMENT LEGEND" ON SHEET A-101.



**STRUCTURAL FLOOR PLAN**

SCALE: 1/4" = 1'-0"



**STRUCTURAL NOTES**

- 1) THE CUSTOMER/CONTRACTOR SHALL FURNISH AND INSTALL ALL STRUCTURAL SUPPORT MEMBERS AND NEEDED HARDWARE FOR THE INSTALLATION OF THE SIEMENS EQUIPMENT.
- 2) THE OVERHEAD STRUCTURAL SUPPORT SYSTEM SHALL BE FIXED, RIGID AND BRACED FOR SWAY.
- 3) ALL STRUCTURAL SUPPORT MEMBERS SHALL BE TRUE, SQUARE, LEVEL, PARALLEL AND COPLANAR WITH RESPECT TO EACH OTHER, WITH A HORIZONTAL STRUCTURAL SUPPORT MEMBER TO BE LOCATED AND SET WITH A TRANSIT.
- 4) ALL STRUCTURAL SUPPORT DETAILS SHOWN ARE SAMPLE DETAILS BASED UPON TYPICAL AND STANDARD BUILDING PRACTICES AND ARE NOT INTENDED AS ACTUAL CONSTRUCTION DETAILS. ALL CONSTRUCTION DETAILS AND SUPPORT CALCULATIONS SHALL BE PREPARED BY A PROFESSIONAL STRUCTURAL ENGINEER AT THE CUSTOMER'S EXPENSE. IN THE EVENT AN EXISTING SUPPORT SYSTEM IS TO BE USED, IT WILL BE THE CUSTOMER'S RESPONSIBILITY TO VERIFY THE INTEGRITY OF THAT SYSTEM.
- 5) MOUNTING PLATES, FRAMES, AND HARDWARE SUPPLIED BY SIEMENS AS DETAILED IN THIS DRAWING SET ARE INSTALLED BY SIEMENS UNLESS OTHERWISE REQUIRED. ANY DEVIATION FROM THE PROVIDED MATERIALS OR MOUNTING METHODS MUST BE DESIGNED AND DOCUMENTED BY THE STRUCTURAL ENGINEER OF RECORD. ALTERNATE MOUNTING MATERIALS (I.E. ANCHORS, THREADED ROD, BACKING PLATES, ETC.) MUST BE SUPPLIED BY THE CUSTOMER/CONTRACTOR. SIEMENS MAY REQUIRE ASSISTANCE FROM THE CUSTOMER/CONTRACTOR WITH INSTALLATION WHEN UTILIZING ALTERNATE MOUNTING MATERIALS.
- 6) ALL CEILING FIXTURES (I.E. AIR SUPPLY GRILLES, AIR RETURN GRILLES, EXHAUST GRILLES, SPRINKLER HEADS, INCANDESCENT AND FLUORESCENT LIGHT FIXTURES, INTERCOM SPEAKERS, MEDICAL GAS COLUMNS, ETC.) SHALL BE INSTALLED FLUSH MOUNTED WITH THE FINISHED CEILING TO PROVIDE FREE AND UNRESTRICTED TRAVEL OF THE SMS CEILING MOUNTED EQUIPMENT.
- 7) THE BOTTOM SIDE OF THE UNISTRUT CEILING GRID AND ANY CEILING MOUNTED SUPPORT PLATES ARE TO BE INSTALLED FLUSH WITH THE FINISHED CEILING. THE CUSTOMER/CONTRACTOR SHALL ALSO PROVIDE COVERSTRIPS FOR THE UNISTRUT.
- 8) THE STRUCTURAL PLANNING AS SHOWN ON THE 1/4" STRUCTURAL PLAN HAS BEEN COORDINATED WITH THE EQUIPMENT LOCATION AS SHOWN ON THE 1/4" EQUIPMENT LAYOUT PLAN. FOR THIS REASON, ANY DEVIATIONS FROM THE STRUCTURAL PLANNING AS SHOWN MUST BE APPROVED BY SMS PLANNING DEPARTMENT.
- 9) THE STRUCTURAL ENGINEER OF RECORD SHALL BE RESPONSIBLE FOR THE DESIGN AND DETAIL OF FLOOR, WALL AND CEILING STRUCTURES IN ACCORDANCE WITH THE STRUCTURAL INFORMATION SHOWN, AND LOCAL GOVERNING BUILDING CODES.
- 10) ALL ANCHORS, SUPPORTS AND BRACES FOR SECURING THE SIEMENS EQUIPMENT ON THE UNDERSIDE OF THE CONCRETE SLAB (WHETHER SUPPLIED BY SIEMENS OR CONTRACTOR) SHALL BE SECURED IN A MANNER TO PREVENT THEM FROM FALLING DURING A DE-INSTALLATION. ALL WORK FOR SECURING THESE MOUNTS SHALL BE BY THE CONTRACTOR.

**FLOOR LOADING**

DESCRIPTION	F STAT MAX	AMPLITUDE FOR F DYN	SUPPORT SURFACE	DIAMETER
F STAT MAX	STATIC FLOOR LOAD DUE TO GANTRY'S OWN WEIGHT			
AMPLITUDE F DYN	DIFFERENCE BETWEEN MINIMUM AND MAXIMUM FLOOR LOADING DURING GANTRY ROTATION			
TABLE OF PARAMETERS				
GANTRY MEASUREMENT POINTS	F STAT MAX (POUNDS)	AMPLITUDE FOR F DYN (POUNDS)	SUPPORT SURFACE	DIAMETER
(A)	598	+/- 67	2 1/2 IN <sup>2</sup>	2 1/8 IN
(B)	710	+/- 67		
(C)	996	+/- 67		
(D)	607	+/- 67		

THE FLOOR STRUCTURE MUST WITHSTAND THE OCCUPIED WEIGHT OF THE GANTRY AND THE INDIVIDUAL CONTACT LOADING AREA. DURING GANTRY INSTALLATION AND LEVELING, THE MAXIMUM POSSIBLE LOAD ON ONE GANTRY FOOT IS 1745 POUNDS (WITH THE GANTRY STANDING ON TWO DIAGONAL FEET).

TOTAL STATIC LOAD, RESULTING IN THE CENTER OF GRAVITY OF THE GANTRY: 2911 POUNDS.

Facility Name: Orlando Health Melbourne Hospital  
 Project Name: Free-standing ER - Viera  
 OPC Project #: 2323980034-19887

**ATTENTION:**

- THIS DRAWING IS DESIGNED TO CONFORM TO FEATURES AND EQUIPMENT REQUIREMENTS PRESENTED AT THE TIME OF THEIR PREPARATION. SINCE BOTH THESE FACTORS ARE SUBJECT TO DESIGN MODIFICATION, THEY ARE NOT TO BE USED FOR CONSTRUCTION PURPOSES.

- IT IS RECOMMENDED THAT THE SIEMENS DRAWINGS BE INCORPORATED WITH THE CONSTRUCTION DOCUMENTS FOR REFERENCE.

- ALL DIMENSIONS SHOWN ON THIS DRAWING ARE FROM FINISHED SURFACES.  
 - THIS DRAWING DOES NOT PROVIDE RADIATION SHIELDING REQUIREMENTS FOR X-RAY AND ASSOCIATED EQUIPMENT. THE CUSTOMER IS RESPONSIBLE FOR CONSULTING WITH A REGISTERED RADIATION PHYSICIST TO SPECIFY RADIATION PROTECTION.

**FINISHED ROOM HEIGHT**

FOR CT GANTRY ONLY	MINIMUM 6'-10 11/16"
FOR CT GANTRY WITH INJECTOR ARM	MINIMUM 7'-6 9/16"
CAREVISION MONITOR/CEILING MOUNT	SEE DETAIL ON S-102 SHEET

SYM	DATE	DESCRIPTION
△	02/04/26	R-1019A VERSION DATED 02/04/26 APPROVED BY CUSTOMER FOR FINALS
-ISSUE BLOCK-		

PROJECT MANAGER: MIKE ROWAN  
 TEL: (813) 431-2205 EXT:  
 FAX: (813) 315-6350  
 EMAIL: MIKE.ROWAN@SIEMENS-HEALTHINEERS.COM

**SIEMENS**

**ORLANDO HEALTH / MELBOURNE FSED**  
 LAKE ANDREW DRIVE, MELBOURNE, FL 32940  
 CT SCAN ROOM - SOMATOM GO-TOP

PROJECT #: **2610370**

THE USE OR REPRODUCTION OF THIS TITLE BLOCK WITHOUT SIEMENS AUTHORIZATION WILL RESULT IN PROSECUTION UNDER FULL EXTENT OF THE LAW.

ALL RIGHTS ARE RESERVED.

SCALE: AS NOTED REF. #:  
 CPQ-1648726/P

DATE: T. KELLEY


SHEET: **S-101**

3 OF 6 DRAWN BY: T. KELLEY

GO-TOP  
REV 30



REFERENCE DOCUMENT - NOT FOR CONSTRUCTION



**ELECTRICAL DIMENSION PLAN**

SCALE: 1/4" = 1'-0"

Facility Name: Orlando Health Melbourne Hospital  
 Project Name: Freestanding ER - Viera  
 OPC Project #: 2323960034-19987

**ATTENTION:**

- THIS DRAWING IS DESIGNED TO CONFORM TO FEATURES AND EQUIPMENT REQUIREMENTS PRESENTED AT THE TIME OF THEIR PREPARATION. SINCE BOTH THESE FACTORS ARE SUBJECT TO DESIGN MODIFICATION, THEY ARE NOT TO BE USED FOR CONSTRUCTION PURPOSES.  
 - THIS SET OF PLANS REPRESENTS A COMPLETE SET OF DETAILS AND SHOULD NOT BE SEPARATED.

- IT IS RECOMMENDED THAT THE SIEMENS DRAWINGS BE INCORPORATED WITH THE CONSTRUCTION DOCUMENTS FOR REFERENCE.



- ALL DIMENSIONS SHOWN ON THIS DRAWING ARE FROM FINISHED SURFACES.  
 - THIS DRAWING DOES NOT PROVIDE RADIATION SHIELDING REQUIREMENTS FOR X-RAY AND ASSOCIATED EQUIPMENT. THE CUSTOMER IS RESPONSIBLE FOR CONSULTING WITH A REGISTERED RADIATION PHYSICIST TO SPECIFY RADIATION PROTECTION.

FINISHED ROOM HEIGHT	
FOR CT GANTRY ONLY	MINIMUM 6'-10 11/16"
FOR CT GANTRY WITH INJECTOR ARM	MINIMUM 7'-6 9/16"
CAREVISION MONITOR/CEILING MOUNT	SEE DETAIL ON S-102 SHEET

PROJECT MANAGER: MIKE ROWAN TEL: (813) 431-2205 EXT: VMAIL: FAX: (813) 315-6350 EMAIL: MIKE.ROWAN@SIEMENS-HEALTHINEERS.COM		<b>SIEMENS</b>	
<b>ORLANDO HEALTH / MELBOURNE FSED</b> LAKE ANDREW DRIVE, MELBOURNE, FL 32940 CT SCAN ROOM - SOMATOM GO-TOP			
THE USE OR REPRODUCTION OF THIS TITLE BLOCK WITHOUT SIEMENS AUTHORIZATION WILL RESULT IN PROSECUTION UNDER FULL EXTENT OF THE LAW.		PROJECT #: <b>2610370</b>	SHEET: <b>E-102</b>
ALL RIGHTS ARE RESERVED.		SHEET 5 OF 6	DRAWN BY: T. KELLEY
SCALE: AS NOTED		REF. #: CPQ-1648726/P	DATE: T. KELLEY

SYM	DATE	DESCRIPTION
△	02/04/26	R-1019A VERSION DATED 02/04/26 APPROVED BY CUSTOMER FOR FINALS
-ISSUE BLOCK-		

GO-TOP  
REV 30



### POWER REQUIREMENTS


SYSTEM	SUPPLY VOLTAGE (VOLTS)	POWER CONSUMPTION (kVA)	SUPPLY IMPEDANCE (m $\Omega$ )	MAIN CIRCUIT BREAKER (AMPS)
GANTRY WITH PATIENT TABLE	3 $\phi$ 480/277Y $\pm$ 10%	SEE BELOW	$\leq$ 300	125

POWER CONSUMPTION OF GANTRY WITH PATIENT TABLE  
 OPERATING FOR 4 SEC.  $\leq$  115 kVA  
 OPERATING FOR 10 SEC.  $\leq$  100 kVA  
 OPERATING FOR 30 SEC.  $\leq$  75 kVA  
 OPERATING FOR 50 SEC.  $\leq$  63 kVA  
 OPERATING FOR 100 SEC.  $\leq$  40 kVA  
 STAND-BY -  $\leq$  3 kVA

IF AN ON-SITE TRANSFORMER IS REQUIRED TO OBTAIN CT OPERATING VOLTAGE, IT MUST HAVE A MINIMUM OUTPUT CAPACITY OF 115kVA AND CHARACTERISTICS TO MAINTAIN SUPPLY VOLTAGE AND IMPEDANCE REQUIREMENTS (TRANSFORMER AND CONDUCTORS).

DO NOT CONNECT ANY EXTERNAL USERS TO THE CT POWER LINES.

THE EXAMINATION ROOM SHOULD BE EQUIPPED WITH AT LEAST ONE EMERGENCY POWER OFF (PANIC) BUTTON.



### CABLE SEPARATION

THIS ELECTRICAL RACEWAY PLAN DEPICTED IN THIS DRAWING IS PLANNED ACCORDING TO SIEMENS SYSTEM REQUIREMENTS AND UL CERTIFICATION OF THIS SYSTEM. ADDITIONAL SEPARATION OF THE SYSTEM CABLE SETS INTO SEPARATE OR PARTITIONED RACEWAYS UNLESS OTHERWISE NOTED IS NOT NECESSARY TO ENSURE SEPARATION OF CIRCUITS. INTERCONNECTING CABLE SETS ARE TESTED AS PART OF THE SYSTEM, AND ARE NOT CONSIDERED PREMISE WIRING.

THE CUSTOMER ASSUMES ALL RESPONSIBILITY AND LIABILITY FOR ANY ADDITIONAL SEPARATION REQUIREMENTS INCLUDING, BUT NOT LIMITED TO: DETERMINING THE NEED FOR ADDITIONAL SEPARATION AND DETERMINING ANY ADDITIONAL ITEMS NEEDED OTHER THAN THOSE IDENTIFIED ON THIS PLAN.

- ### GROUNDING NOTES
- EQUIPMENT GROUNDING CONDUCTOR TO COMPLY WITH THE FOLLOWING:
- 1) SIZE GROUNDING WIRE TO SIEMENS EQUIPMENT PER POWER SCHEDULE REQUIREMENTS.
  - 2) DERIVED FROM THE ELECTRICAL SERVICE, TRANSFORMER OR MAIN DISTRIBUTION PANEL FEEDING THE SIEMENS EQUIPMENT.
  - 3) RUN IN THE SAME CONDUIT, TROUGH OR RACEWAY AS THE PHASE CONDUCTORS.
  - 4) CONTINUOUS, WITH NO BREAKS OR USE OF CONDUIT, CHASSIS OR EARTH AS THE SOLE GROUNDING PATH.
  - 5) BONDED TO CHASSIS AND/OR CONDUIT IN ACCORDANCE WITH THE NEC REQUIREMENTS.
  - 6) MINIMIZE CONNECTIONS OR TERMINALS TO ENSURE CONTINUITY OVER THE LIFE OF THE INSTALLATION.
  - 7) AS A NORM, THERE SHOULD NOT BE ANY CURRENT PRESENCE ON THE GROUND CONDUCTOR, BUT IT IS ACCEPTABLE TO HAVE  $\leq$ 500mA DURING OPERATION OF THE IMAGING EQUIPMENT.

Facility Name: Orlando Health Melbourne Hospital  
 Project Name: Freestanding ER - Viera  
 OPC Project #: 2323980034-1987

# ATTENTION:

- THIS DRAWING IS DESIGNED TO CONFORM TO FEATURES AND EQUIPMENT REQUIREMENTS PRESENTED AT THE TIME OF THEIR PREPARATION. SINCE BOTH THESE FACTORS ARE SUBJECT TO DESIGN MODIFICATION, THEY ARE NOT TO BE USED FOR CONSTRUCTION PURPOSES.

- THIS SET OF PLANS REPRESENTS A COMPLETE SET OF DETAILS AND SHOULD NOT BE SEPARATED.

- IT IS RECOMMENDED THAT THE SIEMENS DRAWINGS BE INCORPORATED WITH THE CONSTRUCTION DOCUMENTS FOR REFERENCE.

- ALL DIMENSIONS SHOWN ON THIS DRAWING ARE FROM FINISHED SURFACES.

- THIS DRAWING DOES NOT PROVIDE RADIATION SHIELDING REQUIREMENTS FOR X-RAY AND ASSOCIATED EQUIPMENT. THE CUSTOMER IS RESPONSIBLE FOR CONSULTING WITH A REGISTERED RADIATION PHYSICIST TO SPECIFY RADIATION PROTECTION.

PROJECT MANAGER: MIKE ROWAN TEL: (813) 431-2205 FAX: (813) 315-6350 EMAIL: MIKE.ROWAN@SIEMENS-HEALTHINEERS.COM			
<b>ORLANDO HEALTH / MELBOURNE FSED</b> LAKE ANDREW DRIVE, MELBOURNE, FL 32940 CT SCAN ROOM - SOMATOM GO-TOP			
THE USE OR REPRODUCTION OF THIS TITLE BLOCK WITHOUT SIEMENS AUTHORIZATION WILL RESULT IN PROSECUTION UNDER FULL EXTENT OF THE LAW. ALL RIGHTS ARE RESERVED.		PROJECT #: <b>2610370</b>	SHEET: <b>E-501</b>
SCALE: AS NOTED REF. #: CPQ-1648726/P		SHEET 6 OF 6 DRAWN BY: T. KELLEY	DATE: T. KELLEY

GO-TOP  
REV 30

REFERENCE DOCUMENT - NOT FOR CONSTRUCTION